

LēClamp™ Ligation Device

Specification

▲ Ligation loop

REF	Open Diam. mm	Height mm
Loop-15B0	15	36
Loop-20B0	20	40
Loop-30B0	30	55
Loop-40B0	40	70

▲ Delivery System

REF	Outer Diam. mm	Effective Length mm	Minimum Working Channel mm
LD-165D0	Φ2.1	1650	2.8
LD-195D0	Φ2.1	1950	
LD-230D0	Φ2.1	2300	

▲ Preassembled

REF	Open Diam. mm	Effective Length mm
Loop15B0-LD165D0	15	1650
Loop20B0-LD165D0	20	1650
Loop30B0-LD165D0	30	1650
Loop40B0-LD165D0	40	1650
Loop15B0-LD195D0	15	1950
Loop20B0-LD195D0	20	1950
Loop30B0-LD195D0	30	1950
Loop40B0-LD195D0	40	1950
Loop15B0-LD230D0	15	2300
Loop20B0-LD230D0	20	2300
Loop30B0-LD230D0	30	2300
Loop40B0-LD230D0	40	2300

▲ Short outer sheath is also available. Please consult with LEOMED representative for details.


LeoMed
 乐奥医疗 为生命健康服务
 Website : www.leomed.cn

Leo Medical Co.,Ltd.

Add : 2 Floor, 10 Building, 18 Huashan Road, Changzhou, China.

Tel : +86 519 8515 0730 Fax: +86 519 8515 0790 E-mail : service@leomed.cn

LēClamp™ Ligation Device



LeoMed
 乐奥医疗 为生命健康服务

GI Endoscopic Accessories ▶